



Clinical Application of Guideline-Complete Liquid Biopsy

Why sequence the tumor genome?

Clinical → to guide therapy

Which drugs to give?

Osimertinib for *EGFR*-mutated
NSCLC


Which drugs not to give?

Cetuximab for *KRAS*-mutated
CRC

Research → disease biology and drug discovery

- Pathway analysis to determine new targets
- Clinical trials testing new targeted agents

Genomic testing is critical for patient care



Genomic testing is critical to getting patients on the correct therapy

- Immunotherapies are not very effective in patients who could be treated with targeted therapies

Traditionally, genomic testing is performed on the tissue biopsy that was used for the cancer diagnosis

Tissue samples are often too small to get patients complete genomic testing

- This can lead to patients being given less effective therapies

Clinical case 1



69-year-old woman with with light smoking history admitted with massive stroke



Evaluation showed lung mass with diffuse nodal and bony metastases and endobronchial ultrasound (EBUS) showed lung adenocarcinoma

Guidelines recommend testing for 9 biomarkers in patients with mNSCLC, prior to 1L treatment


	<i>EGFR</i>	<i>ALK</i>	<i>ROS1</i>	<i>BRAF</i>	<i>RET</i>	<i>MET</i>	<i>NTRK</i>	<i>KRAS</i>	<i>ERBB2</i> (HER2)
Biomarkers with FDA-approved targeted therapy	●	●	●	●	●	●	●		
Emerging biomarkers with promising clinical data								●	●

Current Guideline-Based Recommendations for Genomic Testing: NCCN and CAP/IASLC/AMP^{1,2}

The number of mNSCLC biomarkers with FDA-approved therapies is growing

1. NCCN. NCCN Guidelines Non-Small Cell Lung Cancer. 2018. 2. Lindeman NI et al. J Mol Diagn, 2018.

About 30% of biomarkers in mNSCLC tumors are targetable with a high response rate to targeted therapy



Clinical case 1



69-year-old woman with with light smoking history admitted with massive stroke



Evaluation showed lung mass with diffuse nodal and bony metastases and endobronchial ultrasound (EBUS) showed lung adenocarcinoma




PD-L1 IHC showed 16% expression



Guideline-recommended genomic testing identified *EGFR* exon 19 deletion

For some patients with mNSCLC, 1L immunotherapy is suboptimal

Patients with targetable alterations show higher response rates with targeted therapy vs. immunotherapy



TARGETED THERAPY
50-83% Response rates

IMMUNOTHERAPY
4-26% Response rates


Guidelines and drug labels require genomic results for *EGFR* and *ALK* prior to starting immunotherapy¹

1. National Comprehensive Cancer Network. Non-Small Cell Lung Cancer Guidelines. 2020. 2. Mack PC et al. Cancer 2020. 3. Gettinger S et al. J Clin Oncol. 2016. 4. Peters S et al. J Clin Oncol. 2017. 5. Gainor JF et al. Clin Cancer Res. 2016. 6. Dudnik E et al. J Thorac Oncol. 2018. 7. Guisier F et al. J Thorac Oncol. 2019.

Starting 1L therapy for mNSCLC without complete results has consequences¹

- Toxicities included pneumonitis, colitis, and hepatitis
- All patients with severe irAEs required steroids; nearly all required hospitalization


No severe irAEs were identified when patients received TAGRISSO followed by IO




Clinical case 1



Patient responded to an EGFR TKI for >1 year



Pre-treatment imaging



Post-treatment imaging

Clinical case 2



44-year-old woman presented with persistent cough and weight loss; diagnosed with **stage IV lung adenocarcinoma**



Imaging identified a brain lesion suspected to be metastases



Oncologist ordered SOC tissue testing and liquid biopsy testing at time of first oncology appointment

Clinical case 2



44-year-old woman presented with persistent cough and weight loss; diagnosed with **stage IV lung adenocarcinoma**




Liquid biopsy results reported out in 7 days with *RET* fusion



Patient initiated selpercatinib and showed immediate signs of response; continues to improve 2 months later


Though tissue results also found the *RET* fusion, the results were not returned until 3 weeks after liquid biopsy

Liquid biopsy captures tumor DNA non-invasively



- Tumors shed circulating tumor DNA into the blood through a number of different mechanisms
- A liquid biopsy uses a blood sample to perform genomic testing
- A simple blood draw allows you to test any advanced cancer patient
- Some liquid biopsies can cover relevant cancer genes to help inform treatment decisions
- Receive results faster than tissue

History of cell-free DNA



Liquid biopsy avoids challenges that prevent guideline-complete genomic profiling

TISSUE COMPLEXITIES

LENGTHY PROCESS

- ▶ Results can be unpredictable, may take up to a month or longer, and can be incomplete¹⁻⁵

FINITE RESOURCE

- ▶ Exhausted by histopathology stains and PD-L1 testing^{6,7}

PATIENT BURDEN

- ▶ Repeated invasive tissue biopsies expose patients to potential adverse events^{6,7}

PRACTICE / STAFF BURDEN

- ▶ Significant coordination involving multiple care team members^{6,7}

~1 in 2 patients with NSCLC are unable to get complete genomic profiling results from tissue^{2,8}


Tissue biopsies cause harm & often fail

Invasive biopsies can cause harm

- **1-2%** of lung cancer biopsies result in death¹
- **19%** of lung cancer biopsies result in adverse events²
- Average cost of a lung cancer biopsy to Medicare is **\$14,634** because of complications²

Invasive biopsies can miss the target

- Tissue biopsies often fail^{3,4}



¹ National Lung Screening Trial Research Team 2011 *NEJM*
³ Meric-Bernstam 2015 *Journal of Clinical Oncology*

² Lokhandwala 2016 *Clinical Lung Cancer*
⁴ Sundaresan 2015 *Clinical Cancer Research*



High Concordance (**PPV**) with Tissue

6948 consecutive clinical NSCLC samples

2973 samples with an oncogenic driver mutation detected by liquid biopsy

559 samples with tissue genotyping reports

2430 without tissue genotyping reports



“Validation of a plasma-based comprehensive cancer genotyping assay utilizing orthogonal tissue- and plasma-based methodologies”

181 EGFR L858R

98%

291 EGFR exon 19 del

98%

16 EGFR exon 20 ins

100%

37 ALK/ROS1 fusion

92%

26 KRAS mutations

100%

5 BRAF V600E

100%

3 MET exon 14 skipping

100%

86% response to erlotinib, afatinib, gefitinib

82% response to crizotinib

No further testing needed


63% response to MEK inhibitor

> 50% response to crizotinib


All 3 putative false-positive ALK fusions **responded** to targeted therapy

So PPV for ALK/ROS1 fusions was **actually 100%**


Intratumor Heterogeneity: A Needle or Forceps Biopsy May Not Hit the Right Spot



Spatial and Temporal Tumor Heterogeneity



Temporal Heterogeneity Confounds Plasma to Tissue NGS Concordance: Blinded External Validation Study



Prospective study of 102 Consecutive NSCLC patients

International Association for the Study of Lung Cancer



- ❑ Liquid biopsy is particularly recommended **when tumor tissue is scarce or unavailable, or a significant delay (>14 days) is expected** in obtaining tissue
- ❑ Targeted therapy may reliably be predicated upon blood-based results
- ❑ If a plasma-based “hot-spot” test is negative, repeat testing via comprehensive profiling (such as liquid biopsy) should be pursued; however, **if testing initially performed through comprehensive means is negative, no further testing need be pursued**
- ❑ The cobas test covers only some of the *EGFR* alterations and its sensitivity is outperformed by NGS, which can reach higher levels of sensitivity without diminishing specificity:
“**Multiplex panels using NGS platforms are reliable and preferred** as they detect beyond the common mutations...NGS can reach acceptable levels of sensitivity and optimal levels of specificity”

National Comprehensive Cancer Network



National
Comprehensive
Cancer
Network®

- ❑ “if there is **insufficient tissue** to allow testing for all of *EGFR*, *ALK*, *ROS1* and *BRAF*, repeat biopsy and/or plasma testing should be done”
- ❑ “Plasma-based testing should be considered **at progression on EGFR TKIs** for the T790M mutation”

NILE

“Blood-first” benefit found when comparing liquid biopsy directly to SOC tissue testing¹

282 patients prospectively enrolled



Objective

Compare liquid biopsy to tissue testing for guideline-recommended biomarkers* detection

Clinical follow-up at one year or at disease progression

Study Endpoints

- **Primary:** Detection of guideline-recommended genomic biomarkers
- **Secondary:** Median turn around time (TAT)


*8 Guideline-recommended biomarkers: EGFR, ALK, BRAF, ROS1, RET, ERBB2, and MET (amplification and exon 14 skipping)

1. Leighl et al. Clin Cancer Res. 2019.

NILE

Liquid biopsy as effective as SOC tissue testing¹

Number of patients with an identified guideline-recommended biomarker by testing modality





- Liquid biopsy and tissue testing performed similarly in the detection of guideline-recommended biomarkers (27.3% versus 21.3%; $p < 0.0001$, non-inferiority)
- 98% concordance between liquid biopsy and tissue testing for *EGFR*, *ALK*, *ROS1*, and *BRAF**

* At the time of this study, EGFR, ALK, ROS1 and BRAF were the only biomarkers with FDA-approved targeted therapies for NSCLC
1. Leigh et al. Clin Cancer Res. 2019.

NILE

More mNSCLC patients receive guideline-complete testing with liquid biopsy¹




~3x as many patients tested for *EGFR*, *ALK*, *ROS1*, and *BRAF* mutations* with liquid biopsy (95%) vs. tissue testing (31%)

* At the time of this study, EGFR, ALK, ROS1 and BRAF were the only biomarkers with FDA-approved targeted therapies for NSCLC
1. Leigh et al. Clin Cancer Res. 2019.

NILE

Testing with liquid biopsy identifies more patients with actionable biomarkers - faster¹


20% more patients with biomarkers detected when using **liquid biopsy** first



- More patients with biomarkers identified
 - Due to the limitations of tissue testing, only 31% of patients received complete testing
- Median turn-around-time was significantly faster (9 vs. 15 days; $p < 0.0001$)
 - Liquid biopsy median TAT improved to 7 days by the end of the study

Liquid biopsy + tissue testing nearly doubles number of patients identified with targetable mutations¹

Independent (single site) prospective study of 323 patients comparing liquid biopsy + tissue NGS vs. tissue NGS



44% of patients who were eligible for tissue biopsy were unable to get tissue results due to tissue insufficiency

Testing with Liquid biopsy first improves 1L treatment selection for more mNSCLC patients

“**These results** (Aggarwal et al.), combined with the patient satisfaction with the **relative ease of providing blood** rather than a solid tissue sample, suggest a clinical strategy of **pursuing plasma NGS first**, then tissue NGS if plasma NGS cannot detect relevant mutations.”¹

- Drs. Bishal Gyawali and Jack West (*JAMA Onc* editorial)

Expanded mutation detection at diagnosis





Chart review of 1288 patients
– 68% undergenotyped




Genotyping of ctDNA
increases biomarker
identification by 65%

Liquid for patients
with advanced
solid tumors



Who can benefit from liquid biopsy?



Liquid biopsy identifies actionable alterations across cancers

FDA-approved therapies are available for **solid tumor patients** with *NTRK* fusions or MSI-High status

Lung

EGFR •
ALK •
ROS1 •
BRAF •
MET •
RET •
ERBB2 (HER2)
KRAS

Breast

PIK3CA •
ERBB2 (HER2) •
BRCA1 •
BRCA2 •
NTRK •
MSI •

Colorectal

MSI •
KRAS
NRAS
BRAF •
NTRK •
ERBB2 (HER2)

Prostate

MSI •
BRCA1 •
BRCA2 •
NTRK •

- FDA-approved matched therapy

Clinical case 3



63-year-old man initially diagnosed with Stage 2 **CRC**, referred to a large cancer center after **liver metastases** were detected



No previous genomic testing was recorded for the patient

Oncologist **ordered liquid biopsy** to receive guideline-recommended results within 7 days **while locating archival tissue sample**

Liquid biopsy found a **KRAS exon 4 mutation**

Clinical case 3




The oncologist **excluded anti-*EGFR* therapy** from the patient's treatment plan

The patient was put on **chemo + anti-VEGF therapy** and has had a **strong response for 6+ months**



Pre-treatment imaging



Post-treatment imaging

Clinical case 4



63-year-old woman non-smoker **diagnosed with stage IV breast cancer**
Progressed through several lines of therapy, “out of options”



Patient had ER/PR/HER2 testing performed at diagnosis (several years ago)



No comprehensive genomic testing performed at that time; archival tissue at different institution
Liquid biopsy was ordered and in 7 days, identified *PIK3CA* mutation



The patient was **started on PIQRAY® (alpelisib)** May 2020 and **continues to respond today**

Case-based discussion and Q&A



Patient case template

Clinical History

- *Brief description of patient's clinical history*
- *X*
- *X*

Testing

- *Description of prior genomic testing, if available*
- *Liquid biopsy order and results*
- *X*

Treatment and Response

- *Treatment decision based on liquid biopsy results*
- *Patient outcomes, if available*

Pre-treatment imaging

Scan images if possible

Post-treatment imaging

Scan images if possible

Additional clinical
cases (branded)



Liquid “rescues” a tissue T790M-negative

Initial
Presentation



1st Tissue
Biopsy–Lung
EGFR L858R

Progressed
on Erlotinib After
5 Months



2nd Invasive
Biopsy–Lung
EGFR L858R
*No resistance alterations
found using NGS*

Continued
Progression

Whole brain radiation
and six cycles of
Carboplatin,
Pemetrexed,
Bevacizumab, and
Erlotinib, then
Nivolumab

Guardant360





EGFR L858R
EGFR T790M

Response

Clinical and
Radiographic
Response to 3rd
Generation TKI

Baseline
Pre-Osimertinib





After 4 weeks on
Osimertinib

NSCLC case 1: tissue insufficient for genotyping

CLINICAL CASE

- 58 y/o F non-smoker presents with metastatic NSCLC
- Tissue is insufficient for genotyping despite three biopsy attempts
- Guardant360 demonstrated *EML4-ALK* fusion at 0.06% variant allele fraction
- Crizotinib initiated with significant and durable response




Guardant360 “rescues” a tissue insufficient case

NSCLC case 2: tissue tested for 3 of 7 NCCN genes


CLINICAL CASE

- 77 y/o F non-smoker presents with advanced lung adenocarcinoma
- Tissue *EGFR*, *ALK*, *ROS1* negative, no additional in-house testing available
- Guardant360 ordered
- Carboplatin/Taxol initiated
- Guardant360 identified *BRAF* V600E at 2.1%
- Switched to Trametinib + Dabrafenib
- Near complete response

Guardant360 “rescues” an undergenotyped case



Targeted cancer therapies often work dramatically



But at some point
they fail—what to
do next?

Patient case #2

Clinical History

- 43 y/o female non-smoker diagnosed with stage IV lung adenocarcinoma
- Progressed through chemotherapy and immunotherapy
- Hospice considered

Genomic Testing

- Tissue testing performed at diagnosis found no targetable alterations – FISH negative for *ALK*
- Guardant360 was ordered and in 7 days, identified *EML4-ALK* fusion


Treatment and Response

- Started on crizotinib in April 2018
- Continues to respond; currently no evidence of disease


Appendix (includes branded G360 slides)



Liquid Biopsy – capturing tumor DNA non-invasively



Sources of cell-free DNA (cfDNA)



Key facts about cell-free DNA (cfDNA)

Present in circulation at low concentrations in healthy individuals


1-2 hour half-life in circulation

Higher total cfDNA levels in pregnant women and cancer patients

- Mixture of normal (germline) cfDNA and “other” cfDNA
 - Mean fetal fraction in pregnant women: ~12%
 - Median tumor fraction in cancer patients: ~0.4%

Targeted
“hotspot” tests

NGS optimized for
low fractions of the
target DNA



Uncovering a genomic alteration may not always result in the patient receiving appropriate therapy


- Retrospective study of 5,688 patients with advanced NSCLC documented in the Flatiron Health Database between Jan 1, 2011 and Jul 31, 2016
- All patients received genomic testing
- All patients' received first-line therapy and the treatment decision was examined

What the results found...

<50%

of patients with documented *EGFR* or *ALK* alterations were treated with tyrosine kinase inhibitors in the first line

Severe immune-related adverse events are common with sequential PD-(L)1 blockade and osimertinib




- Identified 126 patients treated with both EGFR TKI and PD-(L)1 blockade.
- 15% of patients treated with sequential PD-(L)1 blockade followed by osimertinib developed a severe irAE.
- Severe irAEs were most common among those who began osimertinib within 3 months
- No irAEs were observed when osimertinib preceded PD-(L)1 blockade or when PD-(L)1 was followed by other EGFR-TKIs.

Real-world issues of undergenotyping in community practices

Results from a review of 814 consecutive NSCLC patients treated at 15 community practices


Genotyping rates for genomic alterations with FDA-approved therapies in NSCLC⁶




- Only 8% of patients were tested for all 7 guideline-recommended alterations (including *MET*, *RET*, and *ERBB2* [HER2])

In prospective studies, Guardant360 had comparable detection of targetable biomarkers vs comprehensive tissue testing


Study 1 Results from 128 patients with advanced NSCLC prior to first-line therapy and at progression¹⁵



Study 2 Results from 185 patients with advanced NSCLC prior to first-line therapy¹⁶




Treat what you find with Guardant360¹⁷




- Objective response was achieved regardless of variant allele fraction

Results were reported significantly faster with Guardant360



- Median turn around time was significantly faster for Guardant360 as compared to tissue testing (9 vs 15 days; $p < 0.0001$)
- Guardant360 median turn around time improved to 7 days over the course of the study

Benefits of Guardant360 at metastatic progression



Do not need to deal with the hassles of tracking down archival tissue

Get results in 7 days

Patients are not subjected to the risks of complications of repeat biopsies


Find treatment or clinical trial recommendations by testing over 70 genes, including MSI-High and *NTRK1* fusions

Avoid the complexities and delays of tissue testing

Tissue workflow for genomic profiling




An unpredictable and complex process that requires communication across multiple stakeholders and **can take >14 days** depending on tissue adequacy



Genomic profiling with Guardant360




A fast and reliable process that delivers **guideline-complete results in 7 days**




Guardant360-detected alterations: patients respond

Pooled response rate of Guardant360 NSCLC studies compared to expected response rate




Considerations when choosing a liquid biopsy



“The effective use of this promising new technology by clinicians hinges upon a shrewd understanding of the test characteristics and validation of a given assay . . .”

The COMPLETE study: Comparison Of Matched Plasma Versus TissuE


- Method comparison
- Real-world patients
- Real-world genotyping
- 7,000 consecutive NSCLC
- Actionable driver mutation focus




*3 of 3 Guardant360+ tissue- ALK fusions were treated with and responded to crizotinib
JCO 2016 34:18_suppl, ASCO Abstract LBA11501, updated JTO, 12(1):S263–S264 for World Congress Lung Cancer Abstract OA06.01 2016

Targeting *BRAF* V600E in NSCLC


Tumor Responses to First-line Dabrafenib + Trametinib in *BRAF* V600E-mutant NSCLC



Poor response to IO in patients with *BRAF* V600E NSCLC




Patients treated based on Guardant360 results responded as expected to targeted therapies



86% of NSCLC patients who were treated with targeted therapy based on Guardant360 results achieved a complete response, partial response, or stable disease

Prevalence of MSI-High*

PREVALENCE OF MSI-HIGH ACROSS CANCER¹



- Guidelines recommend MSI testing for patients with advanced cancers: colorectal, endometrial, gastric and gastroesophageal, pancreatic, and prostate
- Multiple checkpoint inhibitors have FDA approved indications of use across advanced cancers for MSI-High patients

*Please note: MSI status is not reported for specimens originating from New York State or for earlier panel versions.

Guardant360 technical specifications for MSI-High*

Technical Specifications ¹	Mutant Allele Frequency	Sensitivity	Specificity
	≥0.1%	95%	100%

Clinical Validation²: Mutant Allele Frequency: ≥0.4%; Sensitivity: 92%; Specificity: 99%;
Reportable Range: Detected/Not Detected

*Please note: MSI status is not reported for specimens originating from New York State or for earlier panel versions / 1. Analytical Validation / 2. Clinical Validation